



Recombinant Danio rerio Cleft lip and palate transmembrane protein 1-like protein (clptm1l)

Product Code	CSB-CF735313DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6DHU1
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	<p>MFPKTSFTSLIVGVFLLYVLHTCWVMYGIVYTKPCEKRRRAESCISPYLAAKPRL QLSVYT ALRPNADGGHSLIHREEEFDVNTKFEKLVNVSLPKKTRKNGTLYAMVFLHQAG VSPWQDP HQVHLVTQLTTYMLPKPPEISLITGQDEPEKPDQQKQSSDSELDRPVSHWRSR LTLNVVS ENFLFDREALPGDVHRYMRVYQSGKKMIYPLLLFVDELSNRVKDLMEINSSST ELPLTIT YDSIALGKLRFWIHMQDAVYSLQQFGFTEKDADEIKGIFVDTNLYFLALTFFVAA FHLLF DFLAFKNDISFWKHKKSMVGMSSKAVLWRCFSTIVIFLYLLDEQTSLLVLPAGI GSLIE VWKVKKAFKIHVIWRGLTPTFLFGKLDSEKRTEEYDTLAMKYLSYLLYPLCVG GAVYAL VFKYKSWYSWIINSLVNGVYAFGFLFMLPQLFVNYKLKSV AHL PWKAFMYKA FNTFIDD VFAFIITMPTSHRLACFRDDVVFLVYLYQRWLYPVDRSRVNEYGVSYDEKPKG KSHED</p>
Source	in vitro E.coli expression system
Target Names	clptm1l
Protein Names	Recommended name: Cleft lip and palate transmembrane protein 1-like protein Short name= CLPTM1-like protein
Expression Region	1-538
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.