



Recombinant Danio rerio Suppressor of tumorigenicity 7 protein-like (st7I)

Product Code	CSB-CF022768DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A9JRA0
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MADSSGSNPQSPGFTEKLKSWLCWSWTYICALWFAMVLTMVYVLRSPCLKQE TVNAASVF LNTLTPKFYVALTGTSSLISGLILIFEWWYFRKYGTSFIEQVSVSHLRPLLGGVE NSGSA GLFSSVNGDAEPRSNVAECKVWRNPLNLFRGAEYSRYTWWTGKEPLTYDM NLSAQDHQT FFLGDTQQLRPEDSVMQKAWRERNPQARIRAAAYQAIELNRECAAAYVLLAEEE ATTITEA ERLFKQALKSAGKDTNLVYIKRRLAMCARKLGRIKEAVKMMRDLMKEFPLLG MLNIHEN LLEALLELQAYADVQAVLAKYDDISLPKSATICYTSALLKARAVSDKFSPEAASR RGLST AEMNAVEAIHRAVEFNPHVPKYLLEMKSLILPPEHILKRGDSEAVAYAFFHLQH WKRAEG ALNLLHCTWEGTFRIIPYPLEKGHLFYPYPGCTETADRELLPSFHEVSVYPKKE LPFFIL FTAGLCSFCAMLAMLTHQFPELMGVFVKAFFSTLFAPLGGFFADKMESFMPSCL WHQLTNV
Source	in vitro E.coli expression system
Target Names	st7I
Protein Names	Recommended name: Suppressor of tumorigenicity 7 protein-like
Expression Region	1-540
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.