



Recombinant Escherichia coli Putative phosphoethanolamine transferase yhbX (yhbX)

Product Code	CSB-CF337157ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P42640
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	<p>MTVFNKFARTFKSHWLLYLCVIVFGITNLVASSGAHMQRLLFFVLTLVVKRIS SLPLR LLVAAPFVLLTAADMSISLYSWCTFGTTFNDGFAISVLQSDPDEVVKMLGMYIP YLCAFA FLSLLFLAVIIKYDVS LPTKKVTGILLLVISGSLFSACQFAYKDAKNKKAFSPYILA SR FATYTPFFNLNYFALAAKEHQRLLSIANTVPYFQLSVRDTGIDTYVLIVGESVRV DNMSL YGYTRSTTPQVEAQRKQIKLFNQAISGAPYTALSVPLSLTADSVLSHDIHNYPD NIINMA NQAGFQTFWLSSQSAFRQNGTAVTSIAMRAMETVYVRGFDELLPHLSQALQ QNTQQKKL IVLHLNGSHEPACSAYPQSSAVFQPQDDQDACYDNSIHYTDSLLGQVFELLKD RRASVMY FADHGLERDPTKKNVYFHGGREASQQAYHVPMFIWYSPVLGDGVDRTTENNI FSTAYNNY LINAWMGVTKPEQPQTLEEVIAHYKGD SRVVDANHDVFDYVMLRKEFTEDKQ GNPTPEGQ G</p>
Source	in vitro E.coli expression system
Target Names	yhbX
Protein Names	Recommended name: Putative phosphoethanolamine transferase yhbX EC= 2.7.-.-
Expression Region	1-541
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.