



Recombinant Mouse Autophagy-related protein 9A (Atg9a)

Product Code	CSB-CF737885MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q68FE2
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	<p>MAQFDTEYQRLEASYSDSPPGEEDLLVHVAEGSKSPWHHIENLDLFFSRVYNL HQQNGFT CMLIGEMFELMQFLFVVAFTTFLVSCVDYDILFANKMVNHSLHPTEPVKVTLPD AFLPAQ VCSARIQENGLITILVIAGVFWIHRLIKFIYNICCYWEIHSFYLHALRIPMSALPY CTW QEVQARIVQTQKEHQICIKRELTELDIYHRILRFQNYMVALVNKSLPLRFRLP GLGEV VFFTRGLKYNFELILFWGPGSLFLNEWSLKAEYKRGGQRLELAQRLSNRILWIG IANFLL CPLILIQILYAFFSYAEVLKREPGALGARCWSLYGRCYLRFHNELEHELQSRL NRGYKP ASKYMNCFLSPLLTLAKNGAFFAGSILAVLIALTIYDEDVLAVEHVLTTVTLLGV TDTV CRSFIPDQHMVFCPEQLLRVILAHIHYPDHWQGVHLGGVAESHRHTPHSHLL PPPSGPG DHRLLPQLYGRGRGCRHLLLCSDGRSPAWPSSVAVWRADRGLSVPASRGR EDRVVAHAL CHHQSRLAAPS</p>
Source	in vitro E.coli expression system
Target Names	Atg9a
Protein Names	Recommended name: Autophagy-related protein 9A Alternative name(s): APG9-like 1
Expression Region	1-551
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.