



Recombinant Innexin-7 (inx-7)

Product Code	CSB-CF631360CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q21123
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MFVFRVLNTVPYTNRTGAKDLVASIHSFLTSNLLVGLAVLISWKQFGGTPIECM VPLDFT SAWVQYSNNYCWAQPTYFIPFTEELVEQVVDPADVADGITIGNGGNRPRFV KKGGEKIS YYQWMSFFLLFEAACFRLPCFIWKYFASQSGMQVGEILRVASDENNAVPLVKK ANIDALC IHLRGVLRFQKRLKLLKIVPHKILRFLNIKYSAYYVTFIYFAKVAFLNVLQSKLL NK YMLPHDRQQNFGFDMWKTIFYGSTNGNETWRENGVFPRVTLCDFETRDMGN VQMHTVQCV LLLNLFTEKIFVFLWAWYILLTAFTVGNLFSWLFVAVFNETYNEHFILNHLEM CET PFDKD DLKNREHVTRFITLYLGTDGLFLLQLIAQHADVFTTELIAALFKTYIEIEAQRATL KQM NAVLPLLRPNDESQVESGKNTAPSTSHNVRRRGTEQLEKNVKSRRQGSLSLSTQL RPFNSFEE PDQPTKKFDDSSSEDENSKKGSKKPSPTKKKASSKNSPQSSSNSRRPSLAHT ASPAFTHH HEPDSKIPKTAEKKHW
Source	in vitro E.coli expression system
Target Names	inx-7
Protein Names	Recommended name: Innexin-7 Alternative name(s): Protein opu-7
Expression Region	1-556
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.