



Recombinant Innexin-10 (inx-10)

Product Code	CSB-CF632546CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q22549
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MVLA AVL SMLRYVGGSDDRDFVDRLHSYFTCNLLIGLAVLVSFKQFGGKPV ECLVPDIFS SSWEQYAENYCWASDTYYVPTNEPVAGLQSDEKRQRKISYYQWVPFFLLLEA ACFRLPSL LWKYL AGHSGIKINEIVKLSSDPNNIKPDIKRANIKSLTVHLQGALRFHRRLLQKK QIRPH RFLWLFNLPYSAFFVTAMYLCTKFFYLANVCLQLMFMNRFLETDKYKWYGMG ALVDLLNG TTWEQSGMFPRVSLCDFDVRVMGMNMQEHTIQCVLVINIFNEKIFILLWFWYLAL LVFTFG SFFYWLLVSLWRHLNVRFIIRHLEMSDIAFDSSDGAQEKVNRFFISNYLKSDGV FVIRMM TLQSGVIFGTDLVQELWRNFHGSEPQLKRSNSAPKIEEREQWWPAYPSLVNPI NPWRYRD ENQANALRWRRALGANVDNSIATQDLMEKLLPQNAHMRPSDELLENRQPIAV QASYIDDE SDGKLKNEEKQNQQNATQPPYCYTNQNPTPYQNQNQIQNQNQYSNYRTPS LSRGTDSRPVSTATDQTKKQSMSTFK
Source	in vitro E.coli expression system
Target Names	inx-10
Protein Names	Recommended name: Innexin-10 Alternative name(s): Protein opu-10
Expression Region	1-559
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.