



Recombinant Escherichia coli Uncharacterized protein yhjU (yhjU)

Product Code	CSB-CF336329ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P37659
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	<p>MTQFTQNTAMPSSLWQYWRGLSGWNFYFLVKFGLLWAGYLNHFHLLNLVFA AFLLMPLPR YSLHRLRHWIALPIGFALFWHDTWLPGPESIMSQGSQVAGFSTDYLDLVTRFI NWQMIG AIFVLLVAWLFLSQWIRITVFVAILLWLNVLTLAGPSFSLWPAGQPTTTVTTTG GNAAA TVAATGGAPVVGDMPAQTAPPTTANLNAWLNNFYNAEAKRKSTFPSSLPADA QPFELLVI NICSLSWSDIEAAGLMSHPLWSHFDFKFNFSATSYSGPAAIRLLRASCQTS HTNLYQ PANNDCYLFDNLSKLGFTQHLMMGHNGQFGGFLKEVRENGGMQSELMDQT NLPVILLGFD GSPVYDDTAVLNRWLDVTEKDKNRSATFYNTLPLHDGNHYPGVSKTADYKA RAQKFFDE LDAFFTELEKSGRKVMVVVVPEHGGALKGDRMQVSGLRDIPSPSITDVPVGVK FFGMKAP HQGAPIVIEQPSSFLAISDLVVRVLDGKIFTEDNVDWKKLTSGLPQTAPVSENS NAVVIQ YQDKPYVRLNGGDWVPYPQ</p>
Source	in vitro E.coli expression system
Target Names	bcsG
Protein Names	Recommended name: Uncharacterized protein yhjU
Expression Region	1-559
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.