



Recombinant Danio rerio Protein tweety homolog 3 (ttyh3)

Product Code	CSB-CF654883DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q32LT7
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	<p>MAAVVNYSPPWWWNLFHRLPHFNLQFQQTSSDFRPDDSDYQKAVLLLGAAL VCLALDLL FLLFYFWLCCCRRKNHDSFNADCCCTAWCVIIATLVCSAGIAVGFYGNGETC DGVTRLT YSLRHANQTVAGIDKLVSESTSSLNETLMEGLVQLETVYSKQTDYLSIVQKLQG QLDELI NLMLGVFPWSSNDLSLDHLASITEQYDWYRWLGYLGLLLFDVIICLLVLVGLIRN SRSIL IGVCFGLVLTLVISWASLGLFESFAVGASDFCVSPDSYITKVTRENAVINQDILQ YYLKC SMGQTNPFQQKLSGSHKALVEMQDDVSELLRSAIRDFPKTKSNLEEMQGVLN STEVSLLHH LTALVDCRSLHMDYVQALTGLCYDVEGLIYLVLSFVTALMFSSIVCSVPHTW QSKRSE EEDGDETSATLGSRAHDNLYRVHMPSLYSCGSSTYGNEASLPAAHTVSNA PVTEYMSQ NANFQTPRCENTPLIGRESPPPSYTSSMRAKYLATSRPDQPRPTESQNGLEP NMRPDLTS RSAPNSRPNSAIHRPHSAIH</p>
Source	in vitro E.coli expression system
Target Names	ttyh3
Protein Names	Recommended name: Protein tweety homolog 3
Expression Region	1-560
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.