



# Recombinant Inner membrane protein yejM (yejM)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-CF359875EGX  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.  |
| <b>Uniprot No.</b>       | P0AD28   |
| <b>Product Type</b>      | Transmembrane Protein  |
| <b>Immunogen Species</b> | Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)  |
| <b>Sequence</b>          | <p>MVTHRQRYREKVSQMVSWGHWFALFNILLSLVIGSRYLFIADWPTTLAGRIYS<br/> YVSIIGH<br/> FSFLVFATYLLILFPLTFIVGSQRLMRFLSVILATAGMTLLIDSEVFTRFHLHLNP<br/> IVW<br/> QLVINPDENEMARDWQLMFISVPVILLLELVFATWSWQKLRSLTRRRRFARPL<br/> AAFLFIA<br/> FIASHVVYIWADANFYRPITMQRANPLSYPMTARRFLEKHGLLDAQEYQRRLI<br/> EQGNPD<br/> AVSVQYPLSELRYRDMGTGQNVLLITVDGLNYSRFEKQMPALAGFAEQNISFT<br/> RHMSSGN<br/> TTDNGIFGLFYGISPSYMDGILSTRTPAALITALNQQGYQLGLFSSDGFTSPLYR<br/> QALLS<br/> DFSMPVVRTQSDEQTATQWINWLGRYAQEDNRWFSWVSFNGTNIDDSNQQ<br/> AFARKYSRAA<br/> GNVDDQINRVLNLRDSDGKLDNTVVIITAGRGIPLSEEEETFDWSHGHLQVPLV<br/> IHWPGT<br/> PAQRINALTDHTDLMTTLMQRLLHVSTPASEYSQGQDLFNPQRRHYWVTAAD<br/> NDTLAITT PKKTLVLNNGKYRRTYNLRGERVKDEKPQLSLLLQVLTDEKRFIAN</p> |
| <b>Source</b>            | in vitro E.coli expression system  |
| <b>Target Names</b>      | yejM   |
| <b>Protein Names</b>     | Recommended name: Inner membrane protein yejM  |
| <b>Expression Region</b> | 1-586  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |