



Recombinant Human E3 ubiquitin-protein ligase synoviolin (SYVN1)

Product Code	CSB-CF803115HU
Abbreviation	Recombinant Human SYVN1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86TM6
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Human E3 ubiquitin-protein ligase synoviolin (SYVN1)
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MFRTAVMMAASLALTGAVVAHAYYKHKQFYPTVVYLTKSSPSMAVLYIQAFVL VFLLGKVMGKVVFFGQLRAAEMEHELLERSWYAVTETCLAFTVFRDDFSPRFVAL FTLLLFLKCFHWLAEDRVDFMERSPNISWLFHCRIVSLMFLLGILDFLFVSHAYH SILTRGASVQLVFGFEYAILMTMVLTIKIKYVLHSVDLQSENPWDNKAVYMLYTE LFTGFIKVLLYMAFMTIMIKVHTFPLFAIRPMYLAMRQFKKAVTDAIMSRRAIRN MNTLYPDATPEELQAMDNVCIICREEMVTGAKRLPCNHIFHTSCLRSWFQRQQ TCPTCRMDVLRASLPAQSPPPPEPADQGPPPAPHPPPLLPQPPNFPQGLLPP FPPGMFPLWPPMGPFPVPPPPSSGEAVAPPSTSAALSRPSGAATTTAAGT SATAASATASGPGSGSAPEAGPAGFPFPWGMPLPPPFAFPMPVPPA GFAGLTPEELRALEGHERQHLEARLQSLRNIHTLLDAAMLQINQYLTVLASLGP PRPATSVNSTEETATTVAAAASSTSIPSSSEATTPTPGASPPAPEMERPPAPESV GTEEMPEDGEPPDAAELRRRRLQKLESPVAH
Research Area	Cell Biology
Source	in vitro E.coli expression system
Target Names	SYVN1
Expression Region	1-617aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal Twin-Strep-tagged
Mol. Weight	70.8 kDa

**Protein Length**

Full Length

Image

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.