



Recombinant Danio rerio Cleft lip and palate transmembrane protein 1 homolog (clptm1)

Product Code	CSB-CF721864DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6DEL2
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	<p>MAAQETQATQQSSVSNGEVSSNGAAASGQVAQTGATAQDPQQQQQQQPSA WQVIKGVLFRR IFIIWAISSWFRGPPSTPDPNTPAGAPRVPSRNLFPKDTLMDLYVYISQSEIFTD FNTE ELFWYQQDLVYGEWTTGDSGDGCYEYYKELDLSESVKQNGSIYMHYFTKSG FHPDPKRRK GQYRRLSTVHATRMLNKYKRRKFLKTKNLLTGETEADPEMIKRAESHGPPVEIIS HWHPNL TINMVDHTAWVKGSVPPPLDQHVKFDVSGDYPIVYFNDYWNLQKDYYPI NDTLLNLP LRLTYCPLSLWRWQLYAAQSARSPWNFLGEDTYEQSDEDQDSVKVALLETNP YLLGLTIV VSIVHSIFEFLAFKNDIQFWNSRQSLEGLSVRSIIFGVFQSLVLLYILDNETNFV VQVS VFIGLLIDFWKITKMDVRLDRENRIAGIVPRLVFKDKSTYVESSTKIYDDMAFK YLSWL LYPLFGCYAVYSLLYVEHKGWYSWVLSMLYGFLLTFGFITMTPQLFINYKMKS V AHLPPWR MLTYKALNTFIDDLFAFVIKMPMMYRIGCLRDDVFFIYLYQRWIYRVDPNRVN EFGTSG VDHSKDSSAQPAAGETPAAITEKPEGEKKND</p>
Source	in vitro E.coli expression system
Target Names	clptm1
Protein Names	Recommended name: Cleft lip and palate transmembrane protein 1 homolog
Expression Region	1-631
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.