



Recombinant Human Macoilin (TMEM57)

Product Code	CSB-CF836669HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8N5G2
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MKRRNADCSKLRRLPLKRNRITEGIYGSTFLYLKFLVWVWALVLLADFDVLEFRFEY LWPFWL FIRSVYDSFRYQGLAFSVFFVCVAFTSNIICLLFIPIQWLFFAASYVWVQYVWH TERGV CLPTVSLWILFVYIEAAIRFKDLKNFHVLDLCPFAAHCIGYPVVTLGFGFKSYVS YKMRL RKQKEVQKENEFYMQLLQQALPPEQQMLQKQEKEAEEAAKGLPDMDSILIH HNGGIPAN KKLSTTLPEIEYREKGKEKDKDAKKHNLGINNNILQPVDISKIIEIYEMENHINS KRLNN DLVGSTENLLKEDSCTASSKNYKNASGVVNSSPRSHSATNGSIPSSSSKNEKK QKCTSKS PSTHKDLMENCIPNNQLSKPDALVRLEQDIKKLKADLQASRQVEQELRSQISSL SSTERG IRSEMGQLRQENELLQNKLNNAVQMKQKDKQNISQLEKKLKAEQEARSFVEK QLMEEKKR KKLEEATAARAVAFAAASRGECTETLRNRIRELEAEGKKLTMDMKVKEDQIRE LELKVQE LRKYKENEKDTEVLMSALSAMQDKTQHLENSLSAETRIKLDLFSALGDAKRQL EIAQGQI LQKDQEIKDLKQKIAEVMAMPSITYSAATSPLSPVSPHYSSKFVETSPSGLDP NASVYQ PLKK
Source	in vitro E.coli expression system
Target Names	MACO1
Protein Names	Recommended name: Macoilin Alternative name(s): Transmembrane protein 57
Expression Region	1-664
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.