



Recombinant *Saccharomyces cerevisiae* Glycerol uptake/efflux facilitator protein (FPS1)

| | |
|--------------------------|--|
| Product Code | CSB-CF326102SVG |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P23900 |
| Product Type | Transmembrane Protein |
| Immunogen Species | <i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast) |
| Sequence | MSNPQKALNDFLSSSESVHTHDSSRKQSNKQSSDEGRSSSQPSHHHSGGTNN NNNNNNNNNSNNNNNGNDGGNDDDDYDYEMQDYRPSQSARPTPTYVPQY SVESGTAFFIQEVIPSA YINTQDINH KDNGPPSASSNRAFRRPRGQTTVSANVLN IEDFYKNADDAHTIPESHLSRRRSR SRATSNAGHSANTGATNGRTTGAQTNM ESNESPRNVPIMVKPKTLYQNPQTPTVLPSTYHPINKWSSVKNTYLKEFLAEF MGTMVMIIFGSAVVCQVNVAGKIQQDNFNVALDNLNVTGSSAETIDAMKSLTS LVSSVAGGTFDDVALGWAAAVVMGYFCAGGSAISGAHLNPSITLANLVYRGFP LKKVPYYFAGQLIGAFTGALILFIWYKRVLQEAYS DWWWNESVAGMFCVFPKP YLSSGRQFFSEFLCGAMLQAGTFALTDPTYTCLSSDVFP LMMFILIFIINASMAYQ TGTAMNLARDLGPRLALYAVGFDHKMLWVHHHHFFWVPMVGPFGALMGGL VYDVCIYQGHESPVNWSLP VYKEMIMRAWFRPGWKKRNRARRTSDLSDFS YNNDDDEEFGERMALQKTKTKSSISDNENEAGEKKVQFKSVQRGKRTFGGIP TILEEEDSIETASLGATTTDSIGLSDTSSSEDSHYGNAKKVT |
| Source | in vitro E.coli expression system |
| Target Names | FPS1 |
| Protein Names | Recommended name: Glycerol uptake/efflux facilitator protein |
| Expression Region | 1-669 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |