



Recombinant Danio rerio Macoilin-2 (tmem57b)

Product Code	CSB-CF647999DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q2TLY1
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	<p>MKRRNADCSKLRRLPKRNRITEGIYGSTFLYLKFLVWVWALVLLADVFLEFRFEY LWPFWL FIRSVYDSFRYQGLAFSVFFVCVAFTSDIICLLFIPVQWLFFAASTYVWVQYVW HTERGV CLPTVSLWILFVYIEAAIRFKDLKHFHVDLCRPFAAHCIGYPVVTLGFGFKSYVS YKMRL RKQKEVQKENEFYMQLLQQALPPEQQLIQRQEREAEAAAAAAAAAASKSIHDV DSPAQAQ NGSAGGKKPSSNTLPELEYREKERGKNESKKQHNHNQNHHSSTSSSILPSVD NKAQEMEY MENHVNSKRLSSSDLLGSTENLLKDEHSSSSSSSTSSNSNKNYKNASGGGGG GGSSSPRG HGTANGSVPSSSGPSSSASSSSKGRKQKYGGGKNSASHRDPVENCIPNNQ LSKPEALVR LEQDVKKLKADLQASRQTEQDLRSQGLSLGTSERSIRSELGQLRQENELLQNK LHNAVQA KQKDKQTLGQLEKRLKAEQEARAAAEKLLAEKKRKKLEEATAARAVALAAAT RGETES LRRRISELEAECKKLTLDIKVKEDQIRELELKVQELHXYKENEKDTEVLMSALSA MQDKT QHLENSLSAETRIKLDLFSALGDAKRQLEIAQGGQILQKDQEIKDLKQKIAEVMVA MPSVV YSADTGSMTPTVTPHYSSKFMDTSPSGLDPNASVYQPLKK</p>
Source	in vitro E.coli expression system
Target Names	maco1b
Protein Names	Recommended name: Macoilin-2 Alternative name(s): Transmembrane protein 57b
Expression Region	1-699
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.