



Recombinant Alpha-hemolysin translocation ATP-binding protein HlyB (hlyB)

| | |
|--------------------------|---|
| Product Code | CSB-CF318710ENL |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P10089 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Escherichia coli |
| Sequence | <p>MDSCHKIDYGLYALEILAQYHNVSVNPEEIKHRFDTDGTGLGLTSWLLAAKSLE LKVKQV KKTIDRLNFISLPALVWREDGCHFILT KVSKEANRYLIFDLEQRNPRVLEQSEFE ALYQG HIILIASRSSVTGKLAKFDFTWFIPAIIKYRKIFIELVSVFLQLFALITPLFFQVVM D KVLVHRGFSTLNVITVALSVVVVFEIILSGLRITYIFAHSTSRIDVELGAKLFRHLLA LPI SYFESRRVGDVARVRELDQIRNFLTGGALTSVLDLLFSIFFAVMWYYSPKLT LVILFS LPCYAAWSVFISPILRRRLDDKFSRNADNQSFLVESVTAINTIKAMAVSPQMTNI WDKQL AGYVAAGFKVTVLATIGQQGIQLIQKTVMIIINLWLGAHLVISGDLSIGQLIAFNML AGQI VAPVIRLAQIWQDFQQVGISVTRLGDVLNSPTESYHGKLALPEINGDITFRNIRF RYKPD SPVILDNINLSIKQGEVIGIVGRSGSGKSTLTKLIQRFYIPENGQVLIDGHDLALA DPNW LRRQGVVVLQDNVLLNRSIIDNISLANPGMSVEKVIYAAKLAGAHDFISELREGY NTIVG EQGAGLSGGQRQRIAIARALVNNPKILIFDEATSALDYESEHVIMRNMHKICKG RTVIII AHRLSTVKNADRIIVMEKGKIVEQQGKHKELLSEPELSYLYQLQSD</p> |
| Source | in vitro E.coli expression system |
| Target Names | hlyB |
| Protein Names | Recommended name: Alpha-hemolysin translocation ATP-binding protein HlyB |
| Expression Region | 1-707 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.