



Recombinant Mouse Palmitoyltransferase ZDHHC5 (Zdhhc5)

Product Code	CSB-CF845010MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8VDZ4
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MPAESGKRFKPSKYVPVSAAAIFLVGATTLFFAFTCPGLSLNVSPAVPIYNAIMF LFVLA NFSMATFMDPGIFPRAEDEDKEDDFRAPLYKTVEIKGIQVRMKWCATCRFYR PPRCSHC SVCDNCVEEFDHHCPWVNNCIGRRNYRYFFLFLLSLTAHIMGVFGFGLLYVLY HIEELSG VRTAVTMAVMCVAGLFFIPVAGLTGFHVVLVARGRTTNEQVTGKFRGGVNPF TNGCCNNV SRVLCSSPAPRYLGRPKEKTIVIRPPFLRPEVSDGQITVKIMDNGIQGELRRTK SKGSL EITESQSADAEP PPPPKPDLSRYTGLRTHLSLATNEDSSLLGKDSPPTPTMYKY RPGYSS SSTSAAMPHSSSAKLSRGDSLKEPTSIADSSRHPSYRSEPSLEPESFRSPTFG KSFHFD LSSGSRSSSLKSAQGTGFELGQLQSIRSEGTTSTSYKSLANQTRNGSLSYDSL LTPSDSP DFESVQAGPEPDPLGYTSPFLSARLAQQREAERHPRLPTGPPHREPSVR YDNLSRHI VASLQEREKLLRQSPPLAGREEEPGLGDSGIQSTPGSGHAPRTSSSSDDSKR SPLSKTPL GRPAVPRFGKPDGLRSRGLGSPEPGTTAPYLGRSISYSSQKAPSGVSETEEV ALQPLLTP KDEVQLKTTYSKSNGQPKSIGSASPGPGQPPLSSPTRGGVKKVSGVGGTTYEI SV
Source	in vitro E.coli expression system
Target Names	Zdhhc5
Protein Names	Recommended name: Palmitoyltransferase ZDHHC5 EC= 2.3.1.- Alternative name(s): Zinc finger DHHC domain-containing protein 5 Short name= DHHC-5
Expression Region	1-715
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.