




Recombinant Human ATP-binding cassette sub-family D member 2 (ABCD2)

Product Code	CSB-CF866202HU
Abbreviation	Recombinant Human ABCD2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UBJ2
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MTHMLNAAADRVKWSTRSSAAKRAACLVAAYALKTLYPPIIGKRLKQSGHGKKK AAAYPAAENTEILHCTETICEKPSPGVNADFFKQLLELRKILFPKLVTTETGWLC LHSVALISRTFLSIYVAGLDGKIVKSIVEKKPRTFIIKLIKWLMIAIPATFVNSAIRYL ECKLALAFRTRLVDHAYETYFTNQTYKVINMDGRLANPDQSLTEDIMMFSQS VAHLYSNLTKPILDVMLTSYTLIQTATSRGASPIGPTLLAGLVVYATAKVLKACS PKFGKLVAAEEAHRKGYLRYVHSRIIANVEEIAFYRGHKVEMKQLQKSYKALAD QMNLILSKRLWYIMIEQFLMKYVWSSSGLIMVAIPIITATGFADGEDGQKQVMV SERTEAFTTARNLLASGADAIERIMSSYKEVTELAGYARVYNMFVWFDEVKRV GIYKRTAVIQESESHTKNGAKVELPLSDTLAIKGVKIDVDHGIICENVPIITPAGEV VASRLNFKVEEGMHLITGPNCGCKSSLFRILSGLWPVYEGVLYKPPPQHMFY IPQRPYMSLGLSLRDQVIYPDSVDDMHDKGYTDQDLERILHNHLYHIVQREGG WDAVMDWKDVLSGGEKQRMGMARMFYHKPKYALLDECTSAVSIDVEGKIFQ AAKGAGISLLSITHRPSLWKYHHTHLLQFDGEGGWRFEQLDRTAIRLTLSEEKQKL ESQLAGIPKMQQRLNELCKILGEDSVLKTIKNEDETS
Research Area	Cancer
Source	in vitro E.coli expression system
Target Names	ABCD2
Expression Region	1-740aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	89.3 kDa

**Protein Length**

Full Length

Image

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.