



Recombinant Human Probable palmitoyltransferase ZDHHC8 (ZDHHC8)

Product Code	CSB-CF890747HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9ULC8
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	<p>MPRSPGTRLKPAKYIPVATAAALLVGSSTLFFVFTCPWLTRAVSPAVPVYNGIIF LFVLA NFSMATFMDPGVFPRADEDEDKEDDFRAPLYKNVDVRGIQVRMKWCATCHF YRPPRCSHC SVCDNCVEDFDHHCPCWVNNCIGRRNYRYFFLFLLSLSAHMVGVAFLVYVL NHAEGLGA AHTTITMAVMCVAGLFFIPVIGLTGFHVVLVTRGRRTTNEQVTGKFRGGVNPFTF GCCGNV EHVLCSPAPRYVVEPPRLPLAVSLKPPFLRPELLDRAAPLKVKLSDNGLKAGL GRSKSK GSLDRLDEKPLDLGPPLPKIEAGTFSSDLQTPRPGSAESALSVQRTSPPTPA MYKFRPA FPTGPKVPFCGPGEQVPGPDSLTLGDDSIKSLDFVSEPSLDLPDYGPGLHAA YPPSPPL SASDAFSGALRSLKASSRRGGDHVALQPLRSEGGPPTPHRSIFAPHALPNR NGSLSYD SLLNPGSPGGHACPAHPAVGVAGYHSPYLHPGATGDPPRPLPRSFSPVLGPR PREPSPVR YDNLRSRTIMASIQERKDREERERLLRSQADSLFGDSGVYDAPSSYSLQQASVL SEGPRGP ALRYGSRDDLAVAGPGFGGARNPALQTSLSLSSSVSRAPRTSSSSLQADQAS SNAPGPRP SSGSHRSPARQGLPSPGTPHSPSYAGPKAVAFIHTDLPEPPPSLTVQRDHP QLKTPPSK LNGQSPGLARLGPATGPPGPSASPTRHTLVKKVSGVGGTTYEISV</p>
Source	in vitro E.coli expression system
Target Names	ZDHHC8
Protein Names	Recommended name: Probable palmitoyltransferase ZDHHC8 EC= 2.3.1.- Alternative name(s): Zinc finger DHHC domain-containing protein 8 Short name= DHHC-8 Zinc finger protein 378
Expression Region	1-765
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.



Protein Length

full length protein

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.