



Recombinant Geobacillus sp. UPF0059 membrane protein GWCH70_3318 (GWCH70_3318)

Product Code	CSB-CF512850GFJ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C5D9N0
Product Type	Transmembrane Protein
Immunogen Species	Geobacillus sp. (strain WCH70)
Sequence	MKMFIGELVALSMMAFALGMDAFSVALGMGLFRLQLKQIFYIGIMIGLFHIIMPF LGMFL GRFLSYQFGSIASYIGGALLLLGIQMIVTSFKKESDRFVSPMGIGLIFFAFSVSL DSFS VGLSLGIYGVRILLTILLFGFFSTVLTWMGLMLGRHFQQWLGAYSEALGGSILLA FGLKL LFSF
Source	in vitro E.coli expression system
Target Names	mntP
Protein Names	Recommended name: UPF0059 membrane protein GWCH70_3318
Expression Region	1-184
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.