



Recombinant Methanothermobacter thermautotrophicus UPF0059 membrane protein MTH_1812 (MTH_1812)

Product Code	CSB-CF521477MSR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O27840
Product Type	Transmembrane Protein
Immunogen Species	Methanothermobacter thermautotrophicus (strain ATCC 29096 / DSM 1053 / JCM 10044 / NBRC 100330 / Delta H) (Methanobacterium thermoautotrophicum)
Sequence	MDLLSMVLIGVGLAMDAFSISVSRGLALHESETNYALISALSFQGFQAAMPVLG WVSGLE IQRLVSALAPWAAFILLIIGLKMIYESLIMEEEEFIFSYRELLVLSIATSIDAFVGV S FALLDISIWLPVIVIGLITFILSLAGSYIGERVGHIFENRLEALGGLILILIGLKILLEN VSFT
Source	in vitro E.coli expression system
Target Names	mntP
Protein Names	Recommended name: UPF0059 membrane protein MTH_1812
Expression Region	1-184
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.