



Recombinant Protein hdeD (hdeD)

Product Code	CSB-CF360163EOD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AET6
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O157:H7
Sequence	MLYIDKATILKFDLEMLKKHRRAIQFIALLFIVGLLCISFPFVSGDILSTVVGALLI CS GIALIVGLFSNRSHNFWPVLVSGFLVAVAYLLIGYFFIRAPELGIFAIAAFIAGLFCV AGV IRLMSWYRQRSMKGSWLQLVIGVLDIVIAWIFLGATPMVSVTLVSTLVGIELIFS AASLF SFASLFVKQQ
Source	in vitro E.coli expression system
Target Names	hdeD
Protein Names	Recommended name: Protein hdeD
Expression Region	1-190
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.