



Recombinant Escherichia coli O81 Electron transport complex protein RnfA (rnfA)

| | |
|--------------------------|--|
| Product Code | CSB-CF485331EOP |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | B7MVA5 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Escherichia coli O81 (strain ED1a) |
| Sequence | MTDYLLLVFGTVLVNMFVLVKFLGLCPFMGVSKKLETAMGMGLATTFVMTLASI CAWLID TWILIPLNLIYLRTLAFILVIAVVVQFTEMVVRKTSPVLYRLLGIFLPLITTNCAVLG VA LLNINLGHNFLQSALYGFSAAVGFSLVMVLFAAIRERLAVADVPAPFRGNAIALI TAGLM SLAFMGFSGLVKL |
| Source | in vitro E.coli expression system |
| Target Names | rsxA |
| Protein Names | Recommended name: Electron transport complex protein RnfA |
| Expression Region | 1-193 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |