



Recombinant Escherichia coli O139:H28 Cysteine/O-acetylserine efflux protein (eamB)

Product Code	CSB-CF423009EJD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A7ZQ23
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O139:H28 (strain E24377A / ETEC)
Sequence	MTPTLSSAFWTYTLITAMTPGPNNILALSSATSHGFRQSTRVLGMSLGFLIVM LLCAGI SFSLAVIDPAAVHLLSWAGAAAYIVWLAWKIATSPTKEDGLQAKPISFWASFALQ FVNVKI ILYGVTALSTFVLPQTQALSWVVGVSVLLAMIGTFGNVCWALAGHLFQRLFRQ YGRQLNI VLALLLIYCAVRIFY
Source	in vitro E.coli expression system
Target Names	eamB
Protein Names	Recommended name: Cysteine/O-acetylserine efflux protein
Expression Region	1-195
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.