



Recombinant Escherichia coli UPF0056 inner membrane protein yhgN (yhgN)

Product Code	CSB-CF354285ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P67143
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MNEIISA AVL LILIMDPLGNLPIFMSVLKHTEPKRRRAIMVRELLIALLVMLVFLFA GEK ILAFSLRAETVSISGGIILFLIAIKMIFPSASGNSSGLPAGEEPFIVPLAIPLVAGPT I LATLMLLSHQYPNQMGHLVIALLLAWGGTFVILLQSSFLRLLGEKGVNALERL MGLILV MMATQMFLDGIRMWMKG
Source	in vitro E.coli expression system
Target Names	yhgN
Protein Names	Recommended name: UPF0056 inner membrane protein yhgN
Expression Region	1-197
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.