



Recombinant Methanothermobacter thermautotrophicus Uncharacterized protein MTH_215 (MTH_215)

Product Code	CSB-CF519652MSR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O26317
Product Type	Transmembrane Protein
Immunogen Species	Methanothermobacter thermautotrophicus (strain ATCC 29096 / DSM 1053 / JCM 10044 / NBRC 100330 / Delta H) (Methanobacterium thermoautotrophicum)
Sequence	MMYKVRNMRETEVVISDKTANPAPLGLLGFGITTILLNLHNAGLFPINSMILAMG FAYGG IAQILASVMEYRKGNTFGTVAFGSYGLFWWSLVLLLVIPNLKFLETSGTAAASA DPVAMA SYLFMWGLFTLLMFIATLKLKRGIQVIFISLAVLFFLLTAGEITGSALITAVAGYEGI FT GAAAMYVGLAEVINETHGRDILPT
Source	in vitro E.coli expression system
Target Names	MTH_215
Protein Names	Recommended name: Uncharacterized protein MTH_215
Expression Region	1-204
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.