



Recombinant Escherichia coli Leucine efflux protein (leuE)

Product Code	CSB-CF637180EGW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q1RAZ2
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain UTI89 / UPEC)
Sequence	MFAEYGVLYNYWTYLVGAIFIVLVPGPNTLFLVKNVSSGMKGGYLAACGVFIG DAVLMFL AWAGVATLIKTTPIILFNIVRYLGAFYLLYLGSKILYATLKGKNSETKSDEPQYGAI FKRA LILSLTNPKAILFYVSFFVQFIDVNAPHTGISFFILATLLELVSFCYLSFLIISGAFVT Q YIRTKKKLAKVGNLIGLMFVGFARLTLQS
Source	in vitro E.coli expression system
Target Names	leuE
Protein Names	Recommended name: Leucine efflux protein
Expression Region	1-212
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.