



Recombinant Rhopalosiphum padi ATP synthase subunit a (MT-ATP6)

Product Code	CSB-CF858968REAC
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9B6H6
Product Type	Transmembrane Protein
Immunogen Species	Rhopalosiphum padi (Bird cherry-oat aphid) (Aphis padi)
Sequence	MTTNLNFNIFDPSTTIFNLEMNWISTLLIILFMPNFLWILPNRMNWLLFKMFNMLN NEMLM LYKMKKTKSPAFLFISLFMIFILLNFFSLFPYIFSSSSHMVFSVTLAIPFWMFFIIL STC KNTKNMIAHLIPLNTPYIYLAPLMTIIEETMSIIRPMSLSIRLTANLIAGHLLMTLLNFN S LMIIIFIQMFMMIFELCVALIQSYVFSILSSLYSSE
Source	in vitro E.coli expression system
Target Names	MT-ATP6
Protein Names	Recommended name: ATP synthase subunit a Alternative name(s): F-ATPase protein 6
Expression Region	1-217
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.