



Recombinant Oncorhynchus masou ATP synthase subunit a (mt-atp6)

Product Code	CSB-CF656397OEH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q36964
Product Type	Transmembrane Protein
Immunogen Species	Oncorhynchus masou (Cherry salmon) (Masu salmon)
Sequence	FMSPTYLGIPLIAVALTLPWILFPTPSARWLNRLITLQGWFNRFQQLLLPLNL GGHK WAALLTSLMLFLITLNLGLLPYTFTPTTQLSLNMGLAVPLWLATVIIGMRNQPT AALGH LLPEGTPVPLIPVLIIEETISLFIPLALGVRLTANLTAGHLLIQLIATAAFVLLPLMP T VAILTSIVLFLLLTLEIAVAMIQAYVFVLLLSLYLQENV
Source	in vitro E.coli expression system
Target Names	mt-atp6
Protein Names	Recommended name: ATP synthase subunit a Alternative name(s): F-ATPase protein 6
Expression Region	1-219
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.