



Recombinant Bacillus subtilis UPF0700 transmembrane protein yoak (yoak)

Product Code	CSB-CF516475BRJ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O34343
Product Type	Transmembrane Protein
Immunogen Species	Bacillus subtilis (strain 168)
Sequence	MTAAAYRNTLLSLLCLTAGIVDVIGYLSLGHVFTANMTGNIVLLGLAIGKSIQVTV FNSL TALIGFICGVIIATLMVGKAENTLWPSAVTKALALEAFILFVFACLSFYRAFVPVHI LII LMSISMGIQTAAKKLGIAGISSTVLTGTLASLLEDISGRLFFKKQKKTFLRDVLR RALA IILYCVGAIIVALAEPDFYHFIIWVPIVLIFGIMMTAKLKLSGEK
Source	in vitro E.coli expression system
Target Names	yoak
Protein Names	Recommended name: UPF0700 transmembrane protein yoak
Expression Region	1-225
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.