



Recombinant Rhinoceros unicornis ATP synthase subunit a (MT-ATP6)

Product Code	CSB-CF856849RCQ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q96064
Product Type	Transmembrane Protein
Immunogen Species	Rhinoceros unicornis (Greater Indian rhinoceros)
Sequence	MNENLFASFTTPTMMGLPIVILIIMSPSIMFPSPNRLINNRLISIQWLLQLTSKQ MMST HNNKGQWTLMMLSLILFIGSTNLLGLLPHSFTPTTQLSMNLGMAIPLWAGTVL TGFRHK TKASLAHFLPQGTPFLIPMLVIIETISLFIQPVALAVRLTANITAGHLLMHLIGGAT LA LMNISPTTSFITLVLVLLTILEFAVALIQAYVFTLLVSLYLHDNT
Source	in vitro E.coli expression system
Target Names	MT-ATP6
Protein Names	Recommended name: ATP synthase subunit a Alternative name(s): F-ATPase protein 6
Expression Region	1-226
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.