



Recombinant Donkey ATP synthase subunit a (MT-ATP6)

Product Code	CSB-CF015070DK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P92480
Product Type	Transmembrane Protein
Immunogen Species	Equus asinus (Donkey) (Equus africanus asinus)
Sequence	MNENLFASFATPTMMGLPIVILIIMFPSILFPSSNRLINNRLISIQQWLVLQTSKQ MMTI HNNKGQWTLMMLSLILFIGSTNLLGLLPHSFTPTTQLSMNLGMAIPLWAGTVF MGFRHK TKAALAHFLPQGTPIFLIPMLVIIETISLFIQPMALAVRLTANITAGHLLIHLIGGATL A LMDISPSTALITFIILLTILEFAVAMIQAYVFTLLVSLYLHDNT
Source	in vitro E.coli expression system
Target Names	MT-ATP6
Protein Names	Recommended name: ATP synthase subunit a Alternative name(s): F-ATPase protein 6
Expression Region	1-226
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.