



Recombinant Escherichia coli Inner membrane protein yohK (yohK)

Product Code	CSB-CF359870ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AD19
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MMANIWWSLPLTLIVFFAARKLAARYKFLLNPLLVAMVVIIPFLMLTGISYDSYF KGSE VLNDLLQPAVVALAYPLYEQLHQIRARWKSIIITICFIGSVVAMVTGTSVALLMGA SPEIA ASILPKSVTTPIAMAVGGSIGGIPAISAVCVIFVGILGAVFGHTLLNAMRIRTKAAR GLA MGTASHALGTARCAELDYQEGAFSSLALVLCGIITSLIAPFLFPIILAVMG
Source	in vitro E.coli expression system
Target Names	yohK
Protein Names	Recommended name: Inner membrane protein yohK
Expression Region	1-231
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.