



Recombinant Human 3-hydroxyacyl-CoA dehydratase 4 (PTPLAD2)

Product Code	CSB-CF716579HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5VWC8
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MGPLALPAWLQPRYRKNAYLFIYYLIQFCGHSWIFTNMTVRFFSFGKDSMVD TFYAIGLV MRLCQSVSLLELLHIYVGIESNHLLPRFLQLTERIIILFVVITSQEEVQEKYVVCVL FVF WNLLDMVRYTYSMLSVIGISYAVLTWLSQTLWMPYPLCVLAEAFAYQSLPYF ESFGTY STKLPDFLSIYFPYVLKIYLMMLFIGMYFTYSHLYSERRDILGIFPIKKKKM
Source	in vitro E.coli expression system
Target Names	HACD4
Protein Names	Recommended name: 3-hydroxyacyl-CoA dehydratase 4 Short name= HACD4 EC= 4.2.1.- Alternative name(s): Protein tyrosine phosphatase-like protein PTPLAD2 Protein-tyrosine phosphatase-like A domain-containing protein 2
Expression Region	1-232
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.