



# Recombinant Bordetella bronchiseptica Aquaporin Z (aqpZ)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-CF760296BTY  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| <b>Uniprot No.</b>       | Q7WKG2   |
| <b>Product Type</b>      | Transmembrane Protein  |
| <b>Immunogen Species</b> | Bordetella bronchiseptica (strain ATCC BAA-588 / NCTC 13252 / RB50)<br>(Alcaligenes bronchisepticus)   |
| <b>Sequence</b>          | <pre>MQPLFKRCGAEFFGTFWLVLGGCGSAVLAAGVPQVGIGYAGVALAFGLTVLT MAYAVGHISGGHFNPAVTVGLAASGRFGWRDVPYIQAQVVGAIVAAATLASI AQGVAGFDLVASKFAANGYGDHSPGKYSMQAALICEIVLSAGFVFVILGATDK RAPAGFAPIPIGLALTLIHLISIPVTNTSVNPARSTGPALFVGGWALEQLWLFWL APIAGALVGALAYRLVGTTPSAQR</pre>                        |
| <b>Source</b>            | in vitro E.coli expression system  |
| <b>Target Names</b>      | aqpZ   |
| <b>Protein Names</b>     | Recommended name: Aquaporin Z  |
| <b>Expression Region</b> | 1-236  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |