



Recombinant *Oryza sativa* subsp. *japonica* Probable ascorbate-specific transmembrane electron transporter 2 (Os04g0533500, LOC_Os04g45090)

Product Code	CSB-CF800030OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7XMK3
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Sequence	MAAGLGVKAAPFTYVAHALAVAAVMVLVWCISFRGGLAFEADNKNLIFNVHP VLMLIGY IILGSEAIMIYKIFPKLNHDTTKLIHLILHAIAIVLGAVGIYCAFKFHNESGIANLYSL H SWLGIGTISLYGIQWIFGFVAFFYPGAAPHVRRGALPWHVLFGLFVYVLTATA ELGLLE KLTLQSSGLDKYGAEAFNVFTGLVVALFGAAVVVAAPAHVEEPEGYAPIP VN
Source	in vitro E.coli expression system
Target Names	Os04g0533500
Protein Names	Recommended name: Probable ascorbate-specific transmembrane electron transporter 2 EC= 1.-.-. Alternative name(s): Cytochrome b561-2
Expression Region	1-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.