



Recombinant Human Protein lifeguard 4 (TMBIM4)

Product Code	CSB-CF884516HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9HC24
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MADPDPRYPRSSIEDDFNYGSSVASATVHIRMAFLRKVYSILSLQVLLTTVTST VFLYFE SVRTFVHESPALILLFALGSLGLIFALILNRHKYPLNLYLLFGFTLLEALTVAVVVT FYD VYIILQAFILTTTVFFGLTVYTLQSKKDFSKFGAGLFALLWILCLSGFLKFFFYSEI MEL VLAAAGALLFCGFIIYDTHSLMHKLSPEEYVLA AISLYLDIINLFLHLLRFLEAVNK K
Source	in vitro E.coli expression system
Target Names	TMBIM4
Protein Names	Recommended name: Protein lifeguard 4 Alternative name(s): Golgi anti-apoptotic protein Protein S1R Transmembrane BAX inhibitor motif-containing protein 4 Z-protein
Expression Region	1-238
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.