



Recombinant Heme exporter protein C (ccmC)

Product Code	CSB-CF364671SZB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0ABM4
Product Type	Transmembrane Protein
Immunogen Species	Shigella flexneri
Sequence	MWKTLHQLAIPPRLYQICGWFIPLWLAIASVVVLTVGWIWGFAPADYQQGNS YRIIYLH VPAAIWSMGIYASMAVAFIGLVWQMKNANLAVAAMAPIGAVFTFIALVTGSA WGKPMWG TWWVWDARLTSELVLLFLYVGVIALWHAFDDRRLAGRAAGILVLIGVVNLPIIHY SVEWW NTLHQGSTRMQQSIDPAMRSPLRWSIFGFLLSATLTLMRMRNLILLMEKRRP WVSELIL KRGRK
Source	in vitro E.coli expression system
Target Names	ccmC
Protein Names	Recommended name: Heme exporter protein C Alternative name(s): Cytochrome c-type biogenesis protein CcmC
Expression Region	1-245
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.