



Recombinant Drosophila melanogaster Probable Bax inhibitor 1 (CG7188)

Product Code	CSB-CF894268DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9VSH3
Product Type	Transmembrane Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sequence	MADTANYINDRFQTFMNGLGDRYEPYVREHLSKVYMLVGSTAAATAMGAMLQ MRDFLDLG VLAAVATLVVLGLHFYKDDGKNYYTRLGMLYAFGFCSGQTLGPLLYICSINP AIIISA LTGTFVTFISLSLSALLAEQGKYLGLGMLVSVINTMALLSLFNMVFKSYFVQVT QLYVG VFVMAAFIVYDTQNIVEKCRNGNRDVVQHALDLFFDVLSMFRRLLIILTQKEER KQNERR QNKTK
Source	in vitro E.coli expression system
Target Names	BI-1
Protein Names	Recommended name: Probable Bax inhibitor 1 Short name= BI-1
Expression Region	1-245
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.