



Recombinant Escherichia coli UPF0259 membrane protein yciC (yciC)

Product Code	CSB-CF534305ENP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B1ITK0
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
Sequence	MSITAQSVYRDTGNFFRNQFMTILLVSLLCAFITVVLGHVVFSPSDAQLAQLNDG VPVSGS SGLFDLVQNMSPQQILLQASAASTFSGLIGNAILAGGVILIIQLVSAGQRVSA LRAIG ASAPILPKLFILIFLTLLVQIGIMLVVPGIIMAILLALAPVMLVQDKMGIFASMRS SM RLTWANMRLVAPAVLSWLLAKTLLLLFASSFAALTPEIGAVLANTLSNLISAILLIY LFR LYMLIRQ
Source	in vitro E.coli expression system
Target Names	yciC
Protein Names	Recommended name: UPF0259 membrane protein yciC
Expression Region	1-247
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.