



Recombinant UPF0259 membrane protein yciC (yciC)

Product Code	CSB-CF852053EOD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8XCB7
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O157:H7
Sequence	MSITAQSVYRDTGNFFRNQFMTILLVSLLCAFITVVLGHVVFSPSDAQLAQLNDG VPVSGS SGLFDLVQNMSPQQILLQASAASTFSGLIGNAILAGGVILIIQLVSAGQRVSA LRAIG ASAPILPKLFILIFLTLLVQIGIMLVVPGIIMAILLALAPVMLVQDKMGVFASIRS SM RLTWANMRLVAPAVLSWLLAKTLLLLFASSFAALTPEIGAVLANTLSNLISAILLIY LFR LYMLIRQ
Source	in vitro E.coli expression system
Target Names	yciC
Protein Names	Recommended name: UPF0259 membrane protein yciC
Expression Region	1-247
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.