



Recombinant Danio rerio Cytochrome b ascorbate-dependent protein 3 (cybasc3)

Product Code	CSB-CF006324DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A3KPR5
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MRGIVGFYITYLLCLILGIACVVLVVHWNFMYRDGFAWDGSSKNFNWHPVLMV TGMLVLY GNAAVVYRIPLTWGHNKLPWKLLHAGLLLLSFIFSVIGLCAVFNFNHNVHHTANL YSLHSW VGICTAALFTAQWVMGFTSFLLPCTPMAVRAFVKPTHVWMGAMILVLSIVSCIS GINEKL FFVLKETTNGTKPYSALPPEAVAANS LGVIIVAFGLVVLKILSNQMWQRPEPGD DEGVYR PLAYDGS
Source	in vitro E.coli expression system
Target Names	cyb561a3
Protein Names	Recommended name: Cytochrome b ascorbate-dependent protein 3 EC= 1.-.-.-
Expression Region	1-247
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.