



Recombinant Zea mays Aquaporin SIP2-1 (SIP2-1)

Product Code	CSB-CF857689ZAX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9ATM1
Product Type	Transmembrane Protein
Immunogen Species	Zea mays (Maize)
Sequence	MSPAPSRPRIRPWLVVGDALALAAAWVCAGALVKLLVYGGGLGLGGRPEAEAVK VLSLSLVYMFLEAWLEAASGGASYNPLTVLAAALASHGGPAVYLFTAFARIPAQV IGAVLGVKLIQVTFPNVGK GARLSVGAHHGALAEGLATFMVVMVSVTLKKKEM KSFFMKTWITSIWKNTIHLSSDITGGIMNPASAFAWAYARGDHTTFDHLVYW LAPLQATLLGVWAVTFFTKPKKIKEQKVDENKIKKE
Source	in vitro E.coli expression system
Target Names	SIP2-1
Protein Names	Recommended name: Aquaporin SIP2-1 Alternative name(s): Small basic intrinsic protein 2-1 ZmSIP2-1 ZmSIP2;1
Expression Region	1-249
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.