



# Recombinant *Beijerinckia indica* subsp. *indica* ATP synthase subunit a (atpB)

<b>Product Code</b>	CSB-CF015070BSP
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B2IGL1
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Beijerinckia indica</i> subsp. <i>indica</i> (strain ATCC 9039 / DSM 1715 / NCIB 8712)
<b>Sequence</b>	MAEPHEIDPIHQFHIDRIIPLHFLGTDVSFTNAAFFMLVIVALASLVLLAGTRNHS LVPG RLQSIAEVSYEFIASSTLQLSSGRDGMRFPPFVFSIFMFVFLANLIGLVPYTFTVTS QIAV TFGLAMIVIGTVVIYGLIKHGTHFLGIFAPSGVSPLLLPFMIMIEVISFISRPISLSIR L FANMLAGHITLKVMGGFVAGLLGAGSVYALVAPLPLAMVVIFTAFELLVAFLQA YVFTIL TCVYLNDVHPGH
<b>Source</b>	in vitro <i>E.coli</i> expression system
<b>Target Names</b>	atpB
<b>Protein Names</b>	Recommended name: ATP synthase subunit a Alternative name(s): ATP synthase F0 sector subunit a F-ATPase subunit 6
<b>Expression Region</b>	1-253
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.