



Recombinant *Oryza sativa* subsp. *indica* GDT1-like protein 4 (Osl_29993)

Product Code	CSB-CF384919OFF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2YXC7
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Sequence	AAAGDQEDPRGGGDNGTARLDRRTKMFLHAARASDGGATGMEKAGLGLFDA FFASLSMIL VSEIGDETFIIAALMAMRHPKSTVLSGALSALVVMILSTGLGRIVPNLISRKHTN SAAT VLYAFFGLRLLYIAWRSDSKASQKKEIEEVEEKLEAGQGKSTFRRIFSRFCTPIF LESFV LTFLAEWGDRSQIATIALATHKNAVGVAVGATLGHTICTSFAVVGGSMASKIS QGTVAT IGLLFLGFSLSSYFYPPL
Source	in vitro <i>E.coli</i> expression system
Target Names	Osl_29993
Protein Names	Recommended name: GDT1-like protein 4
Expression Region	23-281
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.