



Recombinant Pongo pygmaeus Cytochrome c oxidase subunit 3 (MT-CO3)

Product Code	CSB-CF883045EXP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9T9X6
Product Type	Transmembrane Protein
Immunogen Species	Pongo pygmaeus (Bornean orangutan)
Sequence	MVHQSHAYHMLKPSWPLTGALSALLMTSGLAMWFHSTLLLTGMLTNAL TMYQWWRD VVRESTYQGHHTLPVQKGLRYGMILFITSEVFFFAGFFWAFYHSSLAPTPQLG GHWPPTG ITPLNPLEVPLLNTAVLLASGVSITWAHHSLMENNRTQMIQALLITILLGIYFTLLQ ASE YIEAPFTISDGIYGSTFFMTTGFHGLHVIIGSTFLTVCLSCQLLFHFTSKHHFGFE AAAW YWHFVDVVWFLYVSIYWWGS
Source	in vitro E.coli expression system
Target Names	MT-CO3
Protein Names	Recommended name: Cytochrome c oxidase subunit 3 Alternative name(s): Cytochrome c oxidase polypeptide III
Expression Region	1-261
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.