



Recombinant *Gazella bennettii* Cytochrome c oxidase subunit 3 (MT-CO3)

Product Code	CSB-CF015074GAY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P68531
Product Type	Transmembrane Protein
Immunogen Species	<i>Gazella bennettii</i> (Chinkara) (Indian gazelle)
Sequence	MTHQTHAYHVMNPSWPLTGALSALLMTSGLIMWFHFNSTLLMLGLTTNML TMYQWWRD VIRESTFQGHHTPNVQKGLRYGMILFIISEVLFFTGFFWAFYHSSLAPTPELGG CWPPTG IHPLNPLEVPLLNTSVLLASGVSITWAHHSLMEGNRNHMLQALFITIALGVYFTL LQASE YYEAPFTISDGVYGSTFFVATGFHGLHVIIGSTFLIVCFFRQLKFHFTSSHFGF EAAAW YWHFVDVVWFLFLYVSIYWWGS
Source	in vitro E.coli expression system
Target Names	MT-CO3
Protein Names	Recommended name: Cytochrome c oxidase subunit 3 Alternative name(s): Cytochrome c oxidase polypeptide III
Expression Region	1-261
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.