



# Recombinant *Pelea capreolus* Cytochrome c oxidase subunit 3 (MT-CO3)

<b>Product Code</b>	CSB-CF015074ESW
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O47692
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Pelea capreolus</i> (Gray rhebok)
<b>Sequence</b>	MTHQTHAYHVMNPSWPLTGALSALLMTSGLIMWFHFNSTLLMLGLTTNML TMYQWWRD IIRESTFQGHHTPSVQKGLRYGMILFIISEVLFFTGFFWAFYHSSLAPTPELGGC WPPTG IHPLNPLEVPLLNTSVLLASGVSITWAHHSLMEGNRNHMLQALFITIALGVYFTL LQASE YYEAPFTISDGVYGSTFFVATGFHGLHVIIGSTFLIVCFFRQLKFHFTSNHHFGF EAAAW YWHFVDVVWFLFLYVSIYWWGS
<b>Source</b>	in vitro <i>E.coli</i> expression system
<b>Target Names</b>	MT-CO3
<b>Protein Names</b>	Recommended name: Cytochrome c oxidase subunit 3 Alternative name(s): Cytochrome c oxidase polypeptide III
<b>Expression Region</b>	1-261
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.