



Recombinant Post-GPI attachment to proteins factor 2 (tag-189)

Product Code	CSB-CF630466CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q22141
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MAFGDDDILSIPFKYFVICIGGLPSSALLICVILSLLLHFDQATSTHCEVANWLPSI SAA VSTYTPEKYIWRILIGLHIGPRLVVAIAFRNFFLLGSPLRPLTGHKRLRFLCNLACF LNLL ENFFLLALTSISSSEDHSLHAKCFGGFAICSIIMLLSTWLFNETGRRRTATNLGQ RSHEY KILGAAIFVLCFFLGAYLYWRHNTYCEPGIYTLFALVEYSAVLSNIFFHCTLYYDF HGKN IVLTSSFGGGHYNLLPTQIDKDT
Source	in vitro E.coli expression system
Target Names	tag-189
Protein Names	Recommended name: Post-GPI attachment to proteins factor 2
Expression Region	1-263
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.