



Recombinant Spermidine/putrescine transport system permease protein PotC (potC)

Product Code	CSB-CF360280EGX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AFK7
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Sequence	MIGRLLRGGFMTAIYAYLYIPIIIILIVNSFNSSRFGINWQGFTTKWYSLLMNNDL LQAA QHSLTMAVFSATFATLIGSLTAVALYRYRFRGKPFVSGMLFVMMSPDIVMAIS LLVLFM LLGIQLGFWSLLFSHITFCLPFVVTVYSRLKGFVDMLEAAKDLGASEFTILRKI ILPL AMPAVAAGWVLSFTLSMDDVVSSFTGPSYEILPLKIYSMVKVGVSPEVNAL ATILLVL SLVMVIASQLIARDKTKGNTGDVK
Source	in vitro E.coli expression system
Target Names	potC
Protein Names	Recommended name: Spermidine/putrescine transport system permease protein PotC
Expression Region	1-264
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.