



Recombinant Human TLC domain-containing protein 2 (TLCD2)

Product Code	CSB-CF023584HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6NGC4
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAPTGLLVAGASFLAFRGLHWGLRRLPTPESAARDRWQWWNLCVSLAHSLL SGTGALLGL SLYPQMAADPIHGHPRWALVLVAVSVGYFLADGADLLWNQTLGKTWDLCHH LVVVSCLS TAVLSGHYVGFMSVSLLELNSACLHLRKLKLLSRQAPSLAFSVTSWASLATLA LFRLVP LGWMSLWLFQRHHQVPLALVTLGGIGLVTVGIMSILGIRILVNDVLQSRPHPPS PGHEK TRGTRTRRDNGPVTSSNSTLSLKD
Source	in vitro E.coli expression system
Target Names	TLCD2
Protein Names	Recommended name: TLC domain-containing protein 2
Expression Region	1-264
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.